

Alpha Update™



— \$149 —

The new Alpha Update allows the clinician to create a new bite relationship by building a new lower component. You can choose to change the maxillo-mandibular relationship in any direction: vertical, protrusion or even alter the midline. This gives the clinician freedom to create the best possible result, without limitations. Keep the existing appliance for back-up.


Moses Update is available for any patient of record with any Moses Appliance

Moses Alpha

The Moses

Moses Elite

Here's what the lab needs for an Update:



Upper Impression - *with retainer in place*

- Can seat the retainer on model and take impression or Patient wears retainer and take impression



Bite Registration - *with retainer in place*

- Find the ideal mandibular position, in any position
- Use Moses Bite Shims to insure lip seal and max room for tongue



Lower Impression

- Can just send the original impression / model or Take a new lower impression

Available exclusively from:



Chicago Service Center
773-631-8888

For more information on contraindications, warnings, and possible side effects, please visit our website to download or review additional information.

www.themoses.com

© Moses Appliances LLC 2020